

BeneHeart D30

Defibrillator Monitor

Empowering Every Second





Smart Touch, More Intuitive

Capacitive Full-Touch Screen

With advanced high-definition capacitive touch screen technology and flat UI design, D30 offers optimal experience of clearer display and smoother operation, enabling efficiency at your fingertips.



Optimal visibility

- 8" display with 1024x768 resolution
- Auto-brightness

Confident operation


- Gesture-control touchscreen
- Keep physical knob and buttons for key operations

Wide adaptability

- With the intelligent sensing, you can operate D30 normally even if
- the screen is spilled with liquid
 - wearing up to 5 layers of gloves

Visual AlarmSight™

D30 provides a series of problem-solving support with graphical visualization, not only alerting you to the problems, but also suggesting solutions to help resolve them simply.



Extensive Data Storage

- | | | | | |
|----------------------------|----------------------------------|--------------------------------|---------------------------------------|-----------------------------------|
| 1000 sets
events | 120 hours
ECG waveform | 200 hours
trend data | 1000 sets
auto test reports | 8 hours
voice recording |
|----------------------------|----------------------------------|--------------------------------|---------------------------------------|-----------------------------------|



All in one

360J Manual defib /AED /Pacing /Monitoring: ECG, SPO₂, NIBP, CO₂



Ultra-light

Only 4.2kg(with battery) for exceptional portability



Versatile paddles

External paddles with contact indicator improve usability for clinicians



Three-step operation

Power-Charge-Shock buttons for rapid and reliable operation



Auditory and visual indication

Trace the operation process with the help of sound and light indicators



Build-in recorder

Enable multiple report printouts e.g. waveform, alarm, event, review, test




Rescue Triangle, More Comprehensive

High-quality rescue requires the resuscitation team to maintain the high quality and effectiveness of cardiopulmonary resuscitation during the rescue process. It also relies on timely, regular review and analysis of the rescue data and rescue quality that help identify the highlight for routine rescue training and assessment. Only by combining and advancing the three-in-one solution can the quality of rescue and the survival rate of SCA patients be improved.

Faster Resuscitation

QShock™ -Faster to Shock

In line with the BeneHeart series, D30 is also equipped with new QShock™ technology, taking less than 5s from power on to shock.



Shorter Interruption

Mindray has developed state-of-the-art technologies that can filter out the compression interference for both ALS and BLS, to effectively reduce pauses in recognizing heart rhythm^[1].

Raw ECG with compression interference




Filtered ECG without compression interference



[1] only available when you perform CPR using defib pads or the CPR sensor.

CPR Feedback of both Process & Effect

With the CPR process feedback monitoring the real-time performance of rescuers, and the effect evaluation reflecting patients' responses to the rescue, D30 provides a comprehensive assessment of CPR to obtain a satisfactory resuscitation result.



Structured Debriefing

D30's structured debriefing protocols improve the performance of resuscitation teams in subsequent resuscitation events.

- The key data of each rescue is automatically uploaded to the debriefing system.
- CPR quality and defibrillation data included.
- Not only post-event analysis, but also periodic review.



Hands-on Training

- Training mode in D30 help you get real operation experience.
- Both single and multi-person mode are available.
- Support CPR and defib operation training



Redefine Toughness with Innovation

Preparing for varied scenarios and contingencies, D30 is designed to take on every new challenge in hospital environments anytime, anywhere. D30 empowers clinicians to respond confidently to changing clinical requirements and stay prepared for whatever comes their way.



0.75m drop protection of all six sides



IP55 water-/dust-proof



Bend-resistant cable, withstanding more than 400,000 tests



Resistant to 49 disinfectants



Bedrail hook for convenient nosocomial transhipment

Keep Connected with Greater Efficient


M-Connect™ IT solution integrates the patient data from D30 to avoid manual recording errors and allows flexible viewing of patient data on mobile devices anytime, anywhere. Its simple yet robust network connection is compatible with standard information infrastructure in most hospitals.



M-IoT Device Manager

M-IoT Device Manager obtains comprehensive device data to help biomedical engineers ensure the safety and effectiveness of all equipment at all times.

- Real-time failure monitoring and guidance for timely maintenance, reducing equipment downtime
- Battery status monitoring to ensure patient safety by limiting interruptions
- IP/MAC address for network access control to ensure cybersecurity



Pictures and function descriptions are for reference only. For actual configuration, please consult Mindray's sales staff.